


XX XX
 ID: XX Page: Acq Test Date: May. 30, 2024 DOB: Age: (DOB empty) Report: May. 30, 2024

Acquisition Page





6860 SW 81st Street Miami, Florida 33143
305-668-6102 • 1-800-IHSYSTEMS • www.ihsys.com

XX XX
ID: XX Page: Acq Test Date: May. 30, 2024 DOB: Age: (DOB empty) Report: May. 30, 2024
Acquisition Page

Num	Filename	Int	Ear	Stim.	Type	Swps/Art	Rate	Mode	PP Amp	SNR RN	Gain	Filters
1	NXLA84B.2	84SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.45	0.02 19.407	100	100-1500Hz
2	NXLA84B.3	84SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.42	0.01 19.237	100	100-1500Hz
3	NXLA84B.1	84SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.46	0.01 18.364	100	100-1500Hz
4	NXLA80B.3	80SPL	L	Bone	8000Hz(1)	560/7	11.4	Altr	0.31	0.01 12.462	100	100-1500Hz
5	NXLA80B.1	80SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.38	0.01 11.418	100	100-1500Hz
6	NXLA80B.2	80SPL	L	Bone	8000Hz(1)	560/3	11.4	Altr	0.37	0.01 12.372	100	100-1500Hz
7	NXLA70B.3	70SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.39	0.04 3.307	100	100-1500Hz
8	NXLA70B.2	70SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.57	0.01 3.384	100	100-1500Hz
9	NXLA70B.1	70SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.27	0.01 3.842	100	100-1500Hz
10	NXLA65B.3	65SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.24	0.03 1.643	100	100-1500Hz
11	NXLA65B.2	65SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.56	0.07 2.899	100	100-1500Hz
12	NXLA65B.4	65SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.22	0.04 1.686	100	100-1500Hz
13	NXLA60B.1	60SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.27	0.06 1.686	100	100-1500Hz
14	NXLA60B.3	60SPL	L	Bone	8000Hz(1)	560/2	11.4	Altr	0.68	0.16 1.365	100	100-1500Hz
15	NXLA60B.2	60SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.36	0.11 1.199	100	100-1500Hz
16	NXLA60B.4	60SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.29	0.08 0.894	100	100-1500Hz
17	NXLA55B.2	55SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.46	0.13 0.812	100	100-1500Hz
18	NXLA55B.1	55SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.25	0.11 0.736	100	100-1500Hz
19	NXLA55B.4	55SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.29	0.10 0.779	100	100-1500Hz
20	NXLA55B.3	55SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.23	0.05 0.937	100	100-1500Hz
21	NXLA20B.1	20SPL	L	Bone	8000Hz(1)	560/0	11.4	Altr	0.22	0.46 0.107	100	100-1500Hz
22	NXLA20B.2	20SPL	L	Bone	8000Hz(1)	560/1	11.4	Altr	0.28	0.05 1.375	100	100-1500Hz


System SNR & RN Region: 4.00 - 9.00 ms *-indicates different region used

Num	Int	Ear	Peaks: Latency(ms) Amp(uV) (AR=Amp Ratio)
1	84	L	
2	84	L	
3	84	L	
4	80	L	
5	80	L	
6	80	L	
7	70	L	
8	70	L	
9	70	L	
10	65	L	
11	65	L	
12	65	L	
13	60	L	
14	60	L	
15	60	L	
16	60	L	
17	55	L	
18	55	L	
19	55	L	
20	55	L	
21	20	L	
22	20	L	



SmartEP ASSR OAE VRA
6860 SW 81st Street Miami, Florida 33143
305-668-6102 • 1-800-IHSYSTEMS • www.ihsys.com

XX XX
ID: XX Page: 1 Test Date: May. 30, 2024 DOB: Age: (DOB empty) Report: May. 30, 2024
Page 1



Num Int Ear Peaks: Latency(ms) | Amp(uV) (AR=Amp Ratio)
1 50 L
2 50 L