


XX XX  
ID: XX Page: Acq Test Date: Jun. 03, 2024 DOB: Age: (DOB empty) Report: Jun. 04, 2024  
Acquisition Page







6860 SW 81st Street Miami, Florida 33143  
 305-668-6102 • 1-800-IHSYSTEMS • www.ihsys.com

XX XX  
 ID: XX Page: Acq Test Date: Jun. 03, 2024 DOB: Age: (DOB empty) Report: Jun. 04, 2024  
 Acquisition Page

Num	Filename	Int	Ear	Stim.	Type	Swps/Art	Rate	Mode	PP Amp	SNR RN	Gain	Filters
1	BXLA90B.1	90SPL	L	Bone	500Hz(1)	560/0	11.4	Altr	0.70	0.07 3.487	100	100-1500Hz
2	BXLA90B.3	90SPL	L	Bone	500Hz(1)	560/0	11.4	Altr	0.88	0.09 3.486	100	100-1500Hz
3	BXLA90B.2	90SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.55	0.05 3.657	100	100-1500Hz
4	BXLA85B.3	85SPL	L	Bone	500Hz(1)	560/0	11.4	Altr	0.41	0.08 1.776	100	100-1500Hz
5	BXLA85B.2	85SPL	L	Bone	500Hz(1)	560/0	11.4	Altr	0.69	0.14 2.044	100	100-1500Hz
6	BXLA85B.1	85SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.63	0.10 2.244	100	100-1500Hz
7	BXLA85B.4	85SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.56	0.07 3.125	100	100-1500Hz
8	BXLA80B.2	80SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.95	0.10 5.613	100	100-1500Hz
9	BXLA80B.3	80SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.52	0.10 2.605	100	100-1500Hz
10	BXLA80B.4	80SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.44	0.16 1.149	100	100-1500Hz
11	BXLA80B.5	80SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.53	0.12 1.523	100	100-1500Hz
12	BXLA80B.6	80SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.46	0.12 1.419	100	100-1500Hz
13	BXLA75B.3	75SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.38	0.10 0.993	100	100-1500Hz
14	BXLA75B.2	75SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.45	0.10 0.945	100	100-1500Hz
15	BXLA75B.1	75SPL	L	Bone	500Hz(1)	560/0	11.4	Altr	0.64	0.08 3.571	100	100-1500Hz
16	BXLA20B.2	20SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.41	0.10 1.719	100	100-1500Hz
17	BXLA20B.1	20SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.35	0.06 1.585	100	100-1500Hz

System SNR & RN Region: 4.00 - 9.00 ms \*-indicates different region used

Num	Int	Ear	Peaks: Latency(ms)   Amp(uV) (AR=Amp Ratio)
1	90	L	
2	90	L	
3	90	L	
4	85	L	
5	85	L	
6	85	L	
7	85	L	
8	80	L	
9	80	L	
10	80	L	
11	80	L	
12	80	L	
13	75	L	
14	75	L	
15	75	L	
16	20	L	
17	20	L	



Num	Filename	Int	Ear	Stim.	Type	Swps/Art	Rate	Mode	PP Amp	SNR RN	Gain	Filters
1	BXLA80B.1	80SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	1.17	0.05 10.101	100	100-1500Hz
2	BXLA70B.1	70SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.51	0.20 0.528	100	100-1500Hz
3	BXLA70B.4	70SPL	L	Bone	500Hz(1)	560/1	11.4	Altr	0.41	0.28 0.432	100	100-1500Hz
4	BXLA70B.2	70SPL	L	Bone	500Hz(1)	560/0	11.4	Altr	0.38	0.31 0.490	100	100-1500Hz

System SNR & RN Region: 4.00 - 9.00 ms \*-indicates different region used

Num	Int	Ear	Peaks: Latency(ms)   Amp(uV) (AR=Amp Ratio)
1	80	L	
2	70	L	
3	70	L	
4	70	L	